

therein circumnavigating fractionation effects that lead to complete allele drop-out. Fragment analysis ensues, resulting in an EPG for each cell. Each EPG, therefore, is composed of only a single

Treatment Comparison

Extraction Treatment Variations					
DEParray ¹¹¹ benchtop instrument <u>PP1</u>	Higher Pro K and PBS (1X)				
microelectronics into an image- PP2	Lower Pro K and PBS (1X)				
based system that isolates PP3	Lower Pro K and PBS (0.5X)				
nucleated single-cells.	Lower Pro K and PBS (0.25X)				
Experimental Conditions					
Sample Set	241 Single-cell EPGs(scEPGs)				
Amplification, Separation, AnalysisGlobal Filer™ (29 Cycles), 25 sec injection on 3500 Genetic analyzer, and Osiris (10.3.1) Analysis					

Single-Cell Exploratory Analysis

Exploratory analysis on single leukocytes and epithelial cells collected via two collection techniques: i) semiutomated DEParray[™] NxT and ii) manual pico-pipetting technique.

Ex	perimental Conditions	A. 8000			
Sample Set	Leucocyte : 275 scEPGs Epithelial : 368 scEPGs	5000	Genotypes: 6 No. cells: 275	Genotypes: 5	
Extraction	PicoPure Extraction [Manufacture's recommendation for Epithelial cell and PP3 when collecting with DEPArray [™]]	Unternsity [6000 4000 3000		NO. CEIIS: 300	
Amplification, Separation, Analysis	Global Filer™ (29 Cycles), 25 sec injection on 3500 Genetic analyzer, and Osiris(10.3.1) Analysis	20000 10000 10000	20000 Total 10000		23822







We developed a method to cluster/group profiles by scEPG similarity, wherein the scEPGs of each group are then used to determine the Likelihood Ratio for a cluster of cells. Various methods were evaluated, with the most favorable being the model-based clustering as implemented in the R package mclust [1]. There are two types of clustering errors that might occur: i) overclustering, and ii) misclustering.

Figure (A) is a schematic showing correct clustering, over clustering (known genotype is clustered in more than one group), and misclustering (where at least one group contains more than one genotype).



Conclusion

1. PicoPure[™] DNA extraction treatment coupled with a lower concentration of PBS and Pro K led to higher peak heights, larger allele detection rates, and lower sloping values.

2. MBC, as implemented in *mclust*, showed promising results as it correctly clustered samples consisting of 2- to 5- person mixtures with only one exception.

Reference

1. L. Scrucca, M. Fop, T. B. Murphy, and A. E. Raftery, "Mclust 5: Clustering, classification and density estimation using Gaussian finite mixture models," The R journal, vol. 8, no. 1, pp. 289-317, 2016.

Funding and Collaborations

This work was partially supported by NIJ2020-R2-CX-0032 and NIJ2018-DU-BX-K0185 awarded by the National Institute of Justice, Office of Justice Programs. The opinions, findings, and conclusions or recommendations expressed in this publication are those of the author(s) and do not reflect those of the Department of Justice.

